

Supplementary Table 1. The ratio of FLAIR high-signal intensity and Gd enhancement area (FLAIR/Gd ratio) in terms of the product of the largest diameter and its perpendicular diameter

Case no.	Pre		Post		BRR		Initial
	FLAIR	Gd+	FLAIR	Gd+	FLAIR	Gd+	FLAIR/Gd+ ratio
1	76.32	12.95	74.24	10.56	48.00	3.74	5.89
2	67.50	8.00	33.00	10.08	6.75	4.80	8.44
3	77.52	20.25	44.00	13.30	45.76	6.63	3.83
4	49.80	27.28	34.50	19.50	32.16	20.28	1.83
5	75.26	24.00	47.58	16.38	46.40	14.08	3.14
6	80.60	13.80	63.27	4.50	48.88	15.84	5.84
7	133.44	28.00	116.48	20.28	130.90	21.62	4.77
8	46.98	8.58	12.48	4.20	7.02	5.28	5.48
	37.73	5.50	4.50	1.95	3.12	2.52	6.86
	22.75	4.50	5.00	0.00	5.00	3.80	5.06
9	26.24	3.78	4.68	0.57	8.06	4.14	6.94
10	31.82	2.21	7.77	0.24	7.92	3.20	14.40
	7.04	0.90	0.00	0.20	0.00	0.24	7.82

FLAIR, fluid-attenuated inversion recovery; Gd, gadolinium; BRR, best radiological response